

DePuy/Ace A.R.T. Femoral Screw Insertion Set

Top Tray

Catalog #	Description
1245	Radiographic Ruler
1127	Exchange Tube
14012-14	3.2 mm x 14" Guide Pin (3)
2605	Sleeve Key
2606	Jig Nut
2604	Jig Adaptor
2615	Tommy Bar
2601	Screw Targeting Jig, Distal
1095	Impactor Rod
1096	Sliding Hammer
2616	12.5 mm Cannulated Entry Drill
1291	Guide Wire Grip (2)
2621	Hammer Pad
1098	Solid Extraction Rod

Bottom Tray

Catalog #	Description
1415	Sheath
1416	Trocar
1417	5.3 mm Drill Guide
2608	Slot Targeting Sheath
2629	Short Sheath
2630	3.8 mm Drill Guide, Short
2626	Proximal Targeting Wand
2645	5.3 mm Calibrated Drill Bit, Quick Coupling
14315	Power Adaptor, Quick Coupling
14316	T-Handle, Quick Coupling
1255	6.5 mm Solid Hex Screwdriver Shank
1727	Sliding Hook Depth Gauge
2617	3.8 mm Calibrated 3-Flute Trocar Drill Bit
2619	3.8 mm 3-Flute Trocar Drill Bit
2631	4.5 mm Solid Locking Screwdriver Shank

Storage

Catalog #	Description
2638	Case