

ATS 1000 Cuff Compatibility		
System/Cuff	Cat #	Cuff Description
Single-Use/Quick Connect	DT-7108-01Q	Light blue 9" x 1 ¼" (23 cm x 4.5 cm) single port single bladder
	DT-7112-01Q	Maroon 12" x 2" (30 cm x 5 cm) single port single bladder
	DT-7115-01Q	Gray 15" x 2 ½" (38 cm x 6.5 cm) single port single bladder
	DT-7018-01Q	Orange 18" x 3" (46 cm x 7.5 cm) single port single bladder
	DT-7024-01Q	Yellow 24" x 4" (61 cm x 10 cm) single port single bladder
	DT-7030-01Q	Gray 30" x 4" (76 cm x 10 cm) single port single bladder
	DT-7034-01Q	Dark blue 34" x 4" (86 cm x 10 cm) single port single bladder
	DT-7040-01Q	Green 40" x 4" (102 cm x 10 cm) single port single bladder
	DT-7044-01Q	Green 44" x 4" (112 cm x 10 cm) single port single bladder
	DT-7108-02Q	Light blue 9" x 1 ¼" (23 cm x 4.5 cm) dual port single bladder
	DT-7112-02Q	Maroon 12" x 2" (30 cm x 5 cm) dual port single bladder
	DT-7115-02Q	Gray 15" x 2 ½" (38 cm x 6.5 cm) dual port single bladder
	DT-7218-02Q	Orange 18" x 4" (46 cm x 10 cm) dual port single bladder
	DT-7024-02Q	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder
	DT-7030-02Q	Gray 30" x 4" (76 cm x 10 cm) dual port single bladder
	DT-7034-02Q	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder
	DT-7040-02Q	Green 40" x 4" (102 cm x 10 cm) dual port single bladder
	DT-7044-02Q	Green 44" x 4" (112 cm x 10 cm) dual port double bladder
	DB-7024-02Q	Orange 24" x 5 ½" (61 cm x 14 cm) dual port double bladder
	DT-7108-01L	Light blue 9" x 1 ¼" (23 cm x 4.5 cm) single port single bladder
	DT-7112-01L	Maroon 12" x 2" (30 cm x 5 cm) single port single bladder
	DT-7115-01L	Gray 15" x 2 ½" (38 cm x 6.5 cm) single port single bladder
	DT-7018-01L	Orange 18" x 3" (46 cm x 7.5 cm) single port single bladder
	DT-7024-01L	Yellow 24" x 4" (61 cm x 10 cm) single port single bladder
	DT-7030-01L	Gray 30" x 4" (76 cm x 10 cm) single port single bladder
Single-Use/Luer Lock	DT-7034-01L	Dark blue 34" x 4" (86 cm x 10 cm) single port single bladder
	DT-7040-01L	Green 40" x 4" (102 cm x 10 cm) single port single bladder
	DT-7044-01L	Green 44" x 4" (112 cm x 10 cm) single port single bladder
	DT-7108-02L	Light blue 9" x 1 ¼" (23 cm x 4.5 cm) dual port single bladder
	DT-7112-02L	Maroon 12" x 2" (30 cm x 5 cm) dual port single bladder
	DT-7115-02L	Gray 15" x 2 ½" (38 cm x 6.5 cm) dual port single bladder
	DT-7218-02L	Orange 18" x 3" (46 cm x 7.5 cm) dual port single bladder
	DT-7024-02L	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder
	DT-7030-02L	Gray 30" x 4" (76 cm x 10 cm) dual port single bladder
	DT-7034-02L	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder
	DT-7040-02L	Green 40" x 4" (102 cm x 10 cm) dual port single bladder
	DT-7044-02L	Green 44" x 4" (112 cm x 10 cm) dual port single bladder
	DB-7018-01L	Royal blue 18" x 5 ½" (46 cm x 14 cm) single port double bladder
	DB-7018-02L	Navy blue 18" x 5 ½" (46 cm x 14 cm) dual port double bladder
	DB-7024-01L	Navy blue 24" x 5 ½" (61 cm x 14 cm) single port double bladder
	DB-7024-02L	Navy blue 24" x 5 ½" (61 cm x 14 cm) dual port double bladder
	90-6081Q	Light blue 8" x 2 ½" (20 cm x 6.5 cm) single port single bladder
	90-6121Q	Maroon 12" x 3" (30 cm x 7.5 cm) single port single bladder
	90-6181Q	Orange 18" x 3 ½" (46 cm x 9 cm) single port single bladder
	90-6241Q	Yellow 24" x 4" (61 cm x 10 cm) single port single bladder
	90-6341Q	Dark blue 34" x 4" (86 cm x 10 cm) single port single bladder
	90-6441Q	Green 44" x 4 ½" (102 cm x 11.5 cm) single port single bladder
	90-60181Q	Navy blue 18" x 5 ½" (46 cm x 14 cm) single port double bladder
	90-60241Q	Navy blue 24" x 5 ½" (61 cm x 14 cm) single port double bladder
	90-6082Q	Light blue 8" x 2 ½" (20 cm x 6.5 cm) dual port single bladder

	90-6122Q	Maroon 12" x 3" (30 cm x 7.5 cm) dual port single bladder
	90-6182Q	Orange 18" x 3 ½" (46 cm x 9 cm) dual port single bladder
	90-6242Q	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder
	90-6342Q	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder
	90-6442Q	Green 44" x 4 ½" (112 cm x 11.5 cm) dual port single bladder
	90-60182Q	Navy blue 18" x 5 ½" (46 cm x 14 cm) dual port double bladder
	90-60242Q	Orange 24" x 5 ½" (61 cm x 14 cm) dual port double bladder
Reusable/Luer Lock	90-6081L	Light blue 8" x 2 ½" (20 cm x 6.5 cm) single port single bladder
	90-6121L	Maroon 12" x 3" (30 cm x 7.5 cm) single port single bladder
	90-6151L	Gray 15" x 2 ½" (38 cm x 6.5 cm) single port single bladder
	90-6181L	Orange 18" x 3 ½" (46 cm x 9 cm) single port single bladder
	90-6241L	Yellow 24" x 4" (61 cm x 10 cm) single port single bladder
	90-6341L	Dark blue 34" x 4" (86 cm x 10 cm) single port single bladder
	90-6441L	Green 44" x 4 ½" (102 cm x 11.5 cm) single port single bladder
	90-60181L	Black 18" x 5 ½" (46 cm x 14 cm) single port double bladder
	90-60241L	Navy blue 24" x 5 ½" (61 cm x 14 cm) single port double bladder
	90-6082L	Light blue 8" x 2 ½" (20 cm x 6.5 cm) dual port single bladder
	90-6122L	Maroon 12" x 3" (30 cm x 7.5 cm) dual port single bladder
	90-6152L	Gray 15" x 2 ½" (38 cm x 6.5 cm) dual port single bladder
	90-6182L	Orange 18" x 3 ½" (46 cm x 9 cm) dual port single bladder
	90-6242L	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder
	90-6342L	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder
	90-6442L	Green 44" x 4 ½" (112 cm x 11.5 cm) dual port single bladder
	90-60182L	Navy blue 18" x 5 ½" (46 cm x 14 cm) dual port double bladder
	90-60242L	Orange 24" x 5 ½" (61 cm x 14 cm) dual port double bladder