

Smith & Nephew/Richards AMBI/Classic Compression Hip Screw System

Top Tray Instrumentation

Catalog #	Description
11-0016	Tip Threaded Guide Pin, 3.2 mm (3)
11-0021	Perforation Drill, 6.3 mm
11-0022	AMBI Insertion Wrench (2 pieces)
11-0023	Power Combination Reamer
11-0026	Percutaneous Direct Measuring Gauge
11-0925	Adjustable Angle Guide
11-5036	Guide Pin Placement
11-5037	Fixed Angle Guide, 130°
11-5038	Fixed Angle Guide, 135°
11-5039	Fixed Angle Guide, 140°
11-5041	Fixed Angle Guide, 145°
11-5042	Fixed Angle Guide, 150°
11-5061	Insertion/Removal Wrench (2 pieces)
7111-0014	Lag Screw Tap, must be used with the Quick-Connect T-Handle
7111-0030	Centering Sleeve, for use with Insertion/Removal Wrench
7111-0060	Cannulated Barrel Guide (2 pieces)
7111-0136	Quick-Connect Adaptor – Hall/Jacobs Power
7111-5045	Quick-Connect T-Handle

Bottom Tray Instrumentation

Catalog #	Description
11-0018	Supracondylar Pin Guide, 3.2 mm, 90°
11-0019	Supracondylar Pin Guide, 3.2 mm, 95°
11-0020	Plate Tamper
11-0075	Combination Drill Guide, 3.5 mm
11-0077	Bone Screw Tap, for self-tapping screws
11-0926	Lag Screw Trephine, 9 mm ID
11-5035	Hex Screwdriver
11-5043	Long Slot Drill Guide, 3.5 mm/4.5 mm
11-6500	Obturator
21-0204	Plate Clamp
41-3500	Depth Gauge
7111-0026	Self-Holding Hex Screwdriver
7111-0027	Twist Drill, 4.5 mm
7111-0045	Twist Drill, 3.5 mm
7111-0123	Drill Sleeve, 3.5 mm

Smith & Nephew/Richards AMBI/Classic Compression Hip Screw System

4.5 mm Self-Tapping Cortical Bone Screws

Catalog #	Length	Quantity
12-7026	26 mm	3
12-7028	28 mm	3
12-7030	30 mm	3
12-7032	32 mm	6
12-7034	34 mm	8
12-7036	36 mm	12
12-7038	38 mm	16
12-7040	40 mm	16
12-7042	42 mm	10
12-7044	44 mm	6
12-7046	46 mm	6
12-7048	48 mm	3
12-7050	50 mm	3
12-7052	52 mm	3
12-7054	54 mm	3
12-7056	56 mm	3

7111-5085 Screw Pickup

Storage

Catalog #	Description
7111-0137	4.5 mm Bone Screw Caddy
7111-5090	Case