

Stryker Accolade Basic Instrument Set (1020-9000)

Catalog #	Description
1101-2200	Orthonomic T-Handle
1020-1111	Calcar Planer - Standard
1020-1112	Calcar Planer - Large
1020-1100	Neck Resection Guide
1020-1460	Offset Broach Handle (2)
1601-1210	Modular Box Osteotome
1601-1700	Modular Head Impactor

Head Trials

Catalog #	Description
6264-8-028R	28 mm -4 mm V40 Head Trial, Blue
6264-8-032R	32 mm -4 mm V40 Head Trial, Blue
6264-8-928R	28 mm -2.7mm V40 Head Trial, Light Blue
6264-8-632R	32 mm -2.5 mm V40 Head Trial, Light Purple
6264-8-128R	28 mm+0 V40 Head Trial, Black
6264-8-132R	32 mm+0 V40 Head Trial, Black
6264-8-228R	28 mm+4 V40 Head Trial, Green
6264-8-232R	32 mm +4 mm V40 Head Trial, Green
6264-8-328R	28 mm +8 V40 Head Trial, Gray
6264-8-332R	32 mm +8 V40 Head Trial, Gray
6264-8-428R	28 mm+12, V40 Head Trial, Brown
6264-8-432R	32 mm +12 V40 Head Trial, Brown

Neck Trials

Catalog #	Description
1020-2727	127°, 27 mm Neck Trial
1020-3227	132°, 27 mm Neck Trial
1020-2730	127°, 30 mm Neck Trial
1020-3230	132°, 30 mm Neck Trial
1020-2735	127°, 35 mm Neck Trial
1020-3235	132°, 35 mm Neck Trial
1020-2737	127°, 37 mm Neck Trial
1020-3237	132°, 37 mm Neck Trial
1020-2740	127°, 40 mm Neck Trial
1020-3240	132°, 40 mm Neck Trial

Storage

Catalog #	Description
1020-9000	Single Tray Case
1020-9001	Basic Tray